

Ver.23.04

# One-Step qRT-PCR Kit for probe

(w/o fluorescence dye)

Cat.No. ORQ811-R

**Description:** One-Step qRT-PCR Kit for probe is designed for reverse transcription and quantitative real-time analysis of a specific target RNA by one-step reaction. The kit consists of One-Step RT Enzyme Mix and 2X One-Step Master Mix. The One-Step RT Enzyme Mix contains a thermostable Reverse Transcriptase and a RNase inhibitor. One-Step RT Enzyme Mix can reverse transcribe RNA to cDNA at a wide temperature range from 42 to 60°C and be active against RNase A, RNase B and RNase C. 2X One-Step Master Mix containing specialized Hot Start *Taq* DNA polymerase, which greatly reduce primer-dimer formation and can be activated within 2 minutes, the 2X One-Step Master Mix features high specificity and is suitable for fast cycle program.

## Product contents:

- 1) 2X One-Step Master Mix (TaqMan, no ROX): 500μL
- 2) One-Step RT Enzyme Mix: 100μL

## Storage

- Store at –20°C in dark for 12 months
- Aliquot to avoid multiple freeze-thaw cycles Protect from light

## Recommended primer design

- Amplicon size: 80-150bp
- T<sub>m</sub> value: around 60°C
- Primer length: 17-25mer
- Sequence:
  - 45-55% of GC content is recommended.
  - Avoid regional high GC or AT content.
  - Avoid palindrome sequence.
  - Sequence with G or C at the 3' end is recommended.
- Specificity of primers should be confirmed through a BLAST search.

## Recommended probe design

- T<sub>m</sub> value: 6-10°C higher than primers
- Probe length: 20-30mer
- Sequence:
  - 35-65% of GC content is recommended.
  - Avoid regional high GC or AT content
  - Select the strand contains more C's than G's
  - Avoid palindrome sequence – Avoid a G at the 5' end to prevent quenching of the 5' fluorophore.
- Specificity of probe should be confirmed through a BLAST search.

**Recommended reaction mixture set up for qPCR**

	Volume	Final concentration
Template RNA	Varied	1pg – 1µg
Forward primer (10µM)	Varied	125-900nM
Reverse primer (10µM)	Varied	125-900nM
TaqMan Probe (10µM)	Varied	100-200nM
One-Step RT Enzyme Mix	2µL	1X
2X One-Step Master Mix	10µL	1X
Add D.W to	Up to 20µL	-

**Note:**

- **Template amount varies depending on the copy number of target present in the template solution.**
- **The PCR primer and probe concentration for an optimal qPCR reaction may vary according to primer's and probe's properties. Recommended qPCR program**

**Recommended qPCR Program****Standard program**

Steps	Cycles	Temperature	Time
Reverse transcription	1	42-60°C (42-55°C is recommended)	10 min
Enzyme activation	1	95°C	3 min
Denaturation	40-50	95°C	15 sec
Annealing/Extension		60°C	1 min

**Fast program**

Steps	Cycles	Temperature	Time
Reverse transcription	1	42-60°C (42-55°C is recommended)	5 min
Enzyme activation	1	95°C	20 sec
Denaturation	40-50	95°C	3 sec
Annealing/Extension		60°C	30 sec

**Note**

- You may modify the amount of template, extension time, annealing temperature, and the number of PCR cycles according to the target size, primer's T<sub>m</sub>, and the type of templates for amplification.