

Ver.24021

# h-Taq Polymerase

Cat.No. OHT06

**DESCRIPTION:** h-Taq Polymerase is a thermostable recombinant DNA polymerase derived from thermophilic bacterium *Thermus aquaticus*. Its molecular weight is 94kDa. h-Taq Polymerase can amplify DNA target up to 5kb. The elongation velocity is 2kb/min. It has 5' to 3' polymerase activity but lacks of 3' to 5' exonuclease activity that results in a 3'-dA overhangs PCR product. All components of the 10X h-Taq Buffer are at optimal concentration for efficient amplification. It contributes to highly specific incorporation of primer and template.

## PRODUCT CONTENTS:

- 1) h-Taq Polymerase (5U/ $\mu$ L), 250U
- 2) 10X h-Taq Buffer ( $Mg^{2+}$  Plus), 400 $\mu$ L

## PROTOCOL

1) Add the following components to a sterile microcentrifuge tube sitting on ice:

### a. Recommended PCR assay with 10X h-Taq Buffer ( $Mg^{2+}$ Plus)

Reagent	Volume (50 $\mu$ L)	Final Concentration
10X h-Taq Buffer ( $Mg^{2+}$ Plus)	5 $\mu$ L	1X
dNTPs (10mM each)	1 $\mu$ L	0.2mM each
Forward Primer	- $\mu$ L	0.4-1 $\mu$ M
Reverse Primer	- $\mu$ L	0.4-1 $\mu$ M
h-Taq Polymerase (5U/ $\mu$ L)	0.25-0.5 $\mu$ L	1.25-2.5U/50 $\mu$ L
Template DNA	- $\mu$ L	10pg-1 $\mu$ g
Add D.W to	50 $\mu$ L	

### b. Recommended PCR assay with 10X h-Taq Buffer ( $Mg^{2+}$ free)

Reagent	Volume (50 $\mu$ L)	Final Concentration
10X h-Taq Buffer ( $Mg^{2+}$ free)	5 $\mu$ L	1X
dNTPs (10mM each)	1 $\mu$ L	0.2mM each
Forward Primer	- $\mu$ L	0.4-1 $\mu$ M
Reverse Primer	- $\mu$ L	0.4-1 $\mu$ M
25mM $Mg^{2+}$	- $\mu$ L	1-4mM
h-Taq Polymerase (5U/ $\mu$ L)	0.25-0.5 $\mu$ L	1.25-2.5U/50 $\mu$ L
Template DNA	- $\mu$ L	10pg-1 $\mu$ g
Add D.W to	50 $\mu$ L	

**Table for selection volume of 25mM  $MgCl_2$  solution in a 50 $\mu$ L reaction mix:**

Final $Mg^{2+}$ Con. (mM)	1.0	1.5	2.0	2.5	3.0	4.0
$Mg^{2+}$ (25mM)	2 $\mu$ L	3 $\mu$ L	4 $\mu$ L	5 $\mu$ L	6 $\mu$ L	8 $\mu$ L

**Recommendation amounts of template DNA in a 50µL reaction mix:**

Human genomic DNA	0.1µg – 1µg
Plasmid DNA	0.5ng – 5ng
Phage DNA	0.1ng – 10ng
<i>E.coli</i> genomic DNA	10ng – 100ng

**2) Mix the contents of tube. Cap tubes and centrifuge briefly to collect the contents to the bottom.**

When using a thermal cycler that does not contain a heated lid, overlay the reaction mixture with 25µL mineral oil.

**3) Perform 25-35 cycles of PCR amplification as follows:**

Initial Denaturation	94°C	3 minutes
25-35 Cycles	94°C	30 seconds
	55-68°C	30 seconds
	72°C	1-10 minutes
Final Extension	72°C	10 minutes

**4) Incubate for an additional 10 minutes at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.****5) Analyze the amplification products by agarose gel electrophoresis and visualize by nucleic acid dye staining. Use appropriate molecular weight standards.****APPLICATIONS**

- PCR amplification of DNA fragments as long as 5kb
- DNA labeling
- DNA sequencing
- PCR for cloning

**STORAGE**

Store at -20°C for 12 months.