

Ver.24021

## M-MLV Reverse Transcriptase

Cat.No.ODR23

**Description:** M-MLV Reverse Transcriptase, encoded by Moloney Murine Leukemia Virus is an RNA-dependent DNA polymerase that synthesizes the complementary first strand cDNA from a single-stranded RNA template to which a primer has been hybridized. M-MLV Reverse Transcriptase which is the mutated form of M-MLV Reverse Transcriptase rendering the enzyme thermostable with no RNase H activity. It can retain 100% activity at 50°C; also retains more than 80% activity even at 55°C.

### Product contents:

| Contents                              | Volume |
|---------------------------------------|--------|
| M-MLV Reverse Transcriptase (200U/μL) | 5000U  |
| 5X M-MLV Buffer (50mM DTT included)   | 200μL  |

### Storage

- Store at -20°C

### Features

|                            |          |
|----------------------------|----------|
| RNase H activity           | Negative |
| Concentration              | 200U/μL  |
| Expiration date (at -20°C) | 2 years  |

### Application

The first-strand cDNA synthesis; RT-PCR.

### Note:

Thaw 5X RT Buffer at room temperature just before use and refreeze immediately after use.

### Protocol

First-strand cDNA synthesis (reaction volume 20μL)

- Prepare the following mixture in a PCR tube.

|                     |                                                                                                                                                                |        |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|
| RNA                 | 10ng – 1μg total RNA<br>10ng – 0.5μg mRNA                                                                                                                      | - μg   |
| Primer              | 1μL Oligo (dT) <sub>18</sub> / 1μL Random hexamer<br>or<br>0.5μL Oligo (dT) <sub>18</sub> + 0.5μL Random hexamer<br>or<br>1μL Gene-specific Primer (15 – 20μM) | 1μL    |
| 10mM dNTP           |                                                                                                                                                                | 1μL    |
| RNase-free water to |                                                                                                                                                                | 15.5μL |

**Note:** Prepare aliquots of Oligo(dT)<sub>18</sub> and Random Hexamer to avoid multiple freeze thaw cycles.

2. Heat the mixture at 65°C for 5 minutes and cool down immediately on ice for 1 minute. Spin down the tube by centrifugation briefly.

**Note: Use heating block or PCR machine. DO NOT use water bath.**

3. In the same tube, add the following reaction mixture:

(Total reaction volume- 20µL)

| Reagents                            | Volume |
|-------------------------------------|--------|
| 5X M-MLV Buffer (Included 50mM DTT) | 4µL    |
| M-MLV Reverse Transcriptase         | 0.5µL  |

4. Mix gently by pipetting.

**Note: For random hexamer, incubate the combined reaction mixture at 25°C for 5 minutes and then proceed with the next step. For oligo(dT)<sub>18</sub> or gene-specific primers, directly proceed with the next step.**

5. Incubate at 50°C for 60 minutes.

6. Inactivate the reaction mixture by heating at 95°C for 5 minutes.

**Note: Use heating block or PCR machine. DO NOT use water bath.**

7. The product can be used for PCR immediately or stored at -20°C.

**Note: As a recommendation; cDNA product should compose 10% of total reaction volume of PCR.**