

CD9 Rabbit Polyclonal Antibody

OPRA3419

Host	Rabbit	Storage	-20°C, Avoid freeze / thaw cycles
Isotype	IgG	Applications	WB
Reactivity	N/A	Concentration	0.5mg/mL

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Immunogen Information

Immunogen	Recombinant MRP1(Ser112~Ile195) expressed in <i>E.coli</i> .
Synonyms	MRP1; TSPAN29; 5H9; BA2; BTCC-1; DRAP-27 GIG2; P24; Tetraspanin-29; Cell growth-inhibiting gene 2 protein; Leukocyte antigen MIC3; MIC3

Product Information

Buffer Supplied as solution form in PBS, pH7.4, Containing 0.02%, NaN₃, 50% glycerol.

Applications Western blotting: 0.01-2µg/mL

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures

Legend

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778

Pg-Pig; Z-Zebrafish; X-Xe

