

Ver.260301

# 1X Phosphate Buffered Saline (PBS)

Cat.No.: P1020

## Description

Phosphate Buffered Saline (PBS) is routinely used as a wash buffer in Western blot and immunoprecipitation (IP) procedures, as well as during sample preparation of cell lysates for use in general molecular biology applications.

## Storage

Store at room temperature (15-25°C)

## Applications

- Suitable for use in wash buffers for Western blots and IP, as well as cell lysate sample preparation.
- It has been used for the suspension of bacterial cells.
- It has been used during the fixation and immunostaining procedures

## Notes

- If a precipitate is formed, gently heat the solution till all the precipitate is dissolved.