

MTA2 Rabbit Polyclonal Antibody

OPR19422

Reactivity	H,R	Storage	-20°C, Avoid freeze/ thaw cycles
Host	Rabbit	Applications	WB
Isotype	IgG	Concentration	1 mg/mL

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Immunogen Information

Immunogen	Synthesized peptide derived from human MTA2
Human Swiss-Prot Number	O94776
Synonym	Metastasis-associated 1-like 1 (MTA1-L1 protein) p53 target protein in deacetylase complex

Product Information

Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Dilution	WB: 1:500-2000
Observed MW	73 kDa

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures.

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh- Sheep; Pg-Pig; Z-Zebrafish; X-Xenopus; C-Cow.

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778