

Goat Anti-Rat IgG(H+L) (FITC Conjugated)

OPS1021

Description

| | |
|--------------|--|
| Storage | : -20°C, Avoid freeze / thaw cycles |
| Reactivity | : Rat |
| Immunogen | : Rat IgG |
| Host | : Goat |
| Purification | : Affinity purification |
| Conjugation | : FITC (E _x /E _m :495nm/525nm) |
| Buffer | : PBS with 50% glycerol,pH7.4 |

Applications

IF : 1:20 -100

Research Use

For research use only, not for use in diagnostic procedures.

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh- Sheep; Pg-Pig; Z-Zebrafish;X-Xenopus; C-Cow.

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778