



PROMPT
Proteins for Research



Certificate of Analysis

Product Name : Human Anti-Human-OH Vitamin D Monoclonal Antibody
Catalog Number : 25-OHVD101
Lot Number: 20230613
Source: Synthetic
Form: Liquid
Purity : $\geq 90\%$ (by SDS-PAGE)
Concentration : 5.2 mg/mL
Buffer : 0.01 M PBS (pH 7.4)
Preservative: 0.05% Proclin- 300
Appearance: Clear or slightly opalescent solution
Application: ELISA,CLIA
Storage: Store at $-20 \pm 5^{\circ}\text{C}$ for two years.
Avoid repeated freeze- thaw

Note : For diagnostic use only.Not intended
For therapeutic use.

Quality Department

Date: 2024.06.13



For further details, contact us at

Origin Diagnostics and Research, XIII/712, Origin Building, Clappana, Karunagappally, Kollam - 690525
info@originlab.in | +91 7736957333 | www.originlab.in