

Goat Anti-Chicken IgG - HRP

SE135

Host	Goat	Storage	-20°C, Avoid freeze / thaw cycles
Isotype	IgG	Concentration	1mg/mL

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Buffer	Commercial stabilizer
Dilution	1:500 – 1:20000
Purity	G Protein

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778