

Rabbit Anti-Goat IgG -HRP

OPS004

Host	Rabbit	Storage	-20°C, Avoid freeze / thaw cycles
Isotype	IgG	Applications	WB, IHC, ELISA

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Buffer	0.01M PBS (pH 7.4), with 0.03% Proclin300, and 50% glycerol and 1% BSA.
Dilution	WB 1:1000-10000, IHC 1:50-500, ELISA 1:1000-10000 Use TMB as a substrate (OD450nm) P/N=2.5~3.0 for positive.

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Important Note

1. The coated antibody or antigen is uncomfortable, and it is easy to produce small positive values, or produce mutations, that is, tube skipping.
2. Coating with antiserum that is not purified or the degree of purification is too low, the positive value is small, or there is no positive gradient.
3. Poor sealing or over-concentration of markers can produce higher control values.

Background

HRP conjugated antibody is a kind of immunochemistry reagents widely used in immunology, molecular biology and various branches of clinical medicine. The company prepared a variety of HRP conjugated antibody based on modified sodium periodate oxidation method. 1mL product includes about 1.5-2.5mg HRP and 3-5mg IgG.

Research Use

For research use only, not for use in diagnostic procedures.

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh-Sheep; Pg-Pig; Z-Zebrafish; X-Xenopus; C-Cow.

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778