

Goat Anti-Rabbit IgG -HRP

OPS003

Host	Goat	Storage	-20°C, Avoid freeze / thaw cycles
Isotype	IgG	Applications	WB, IHC, ELISA

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Buffer	PBS with 0.02% sodium azide, and 50% glycerol pH 7.4.
Dilution	WB 1:5000-10000, IHC 1:100-1:200, ELISA: Use at an assay dependent concentration

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures.

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh-Sheep; Pg-Pig; Z-Zebrafish; X-Xenopus; C-Cow.

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778