Goat Anti-Mouse IgG - HRP

OPS002

Host	Goat	**Storage**	-20°C, Avoid freeze / thaw cycles
Isotype	IgG	**Applications**	WB, IHC, ELISA

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Dilution	WB 1:5000-1:10000, IHC 1:100-1:200, ELISA: Use at an assay dependent concentration

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh- Sheep; Pg-Pig; Z-Zebrafish; X-Xenopus; C-Cow.